

## FCPC C - ULTRA FAST CPC



### NEW EQUIPMENT

25 ml or 50 ml column capacity

Fast separation

Minimum sample quantities required



DATA	25ml	50ml
SAMPLE QUANTITY	10 - 500 mg	10 - 1000 mg
FLOW RATE	1- 10 ml/min	
SEPARATION TIME	5 – 30 min	
PRESSURE MAX.	80 bars / 1160 psi	
ROTATION SPEED MAX.	3000 rpm	
SIZE / WEIGHT	300 x 300 x h 555 mm / 40 kg	
POWER	750W – 110 or 220 VAC	
MATERIAL OF CONSTRUCTION	Rotor: stainless steel + PTFE. Wetted parts and frame : stainless steel	

Kromaton equipments are suitable for industrial/laboratory service, with reliable and robust operation, and are compatible with current European standards.